




ProLipo *PLUS*™

Less FAT. More PATIENTS. | High-Power, Laser-Assisted Lipolysis



 SCITON®



Less FAT. More PATIENTS.

ProLipo *PLUS*™ is a minimally invasive, ultra high-powered laser-assisted lipolysis system for melting unwanted fat and coagulating tissue resulting in a tighter, slimmer body.

ProLipo *PLUS* employs dual wavelength capability (1064 nm and 1319 nm), providing the choice to use each wavelength individually or to combine the benefits of both for optimal body sculpting results.

With up to 40 watts of power on 1064 nm and up to 40 watts on 1319 nm, the system provides flexibility in tailoring advanced treatments.

A wide range of customizable treatment parameters enables body sculpting treatments that are personalized for each patient.

ProLipo *PLUS* brings you:

New Revenue Opportunity

Laser-assisted lipolysis is one of today's most popular procedures

Superior Results

Deliver speed and results with best-in-class power output and dual wavelength capability

Unprecedented Flexibility

Personalize treatments with a complete range of precisely adjustable settings

Maximum ROI

Protect your investment with a fully expandable, modular platform and non-proprietary fibers

"I was able to pay for the ProLipo PLUS laser in the first four months of ownership. And, with no recurring disposable cost to deal with, I am enjoying a significant return on my investment."

— Marc Salzman, MD, FACS

"You will be hard-pressed to attain the same degree of skin tightening with more conventional forms of liposuction."

— Neil C. Goodman, MD

Neck



Before



After 6 weeks

Photos courtesy of Roger S. Hogue, MD, RVT

Abdomen



Before



After 2 months

Waist and Thighs



Before



After 3 months

Photos courtesy of Brian V. Heil, MD, FACS

Arm



Before



After

Photos courtesy of Kevin Tehrani, MD

CHIN • JOWLS • NECK • MALE CHEST • ARMS • BACK • ABDOMEN • FLANKS • HIPS • WAIST • BUTTOCKS • THIGHS • KNEES

ProLipo *PLUS* brings your patients:

Immediately Visible Results

Eliminates fat and tightens skin through heating and coagulation of soft tissue with results that the patient can see immediately.

Comfortable and Quick Body Contouring

ProLipo treatment with ProLipo *PLUS* is minimally invasive, a quick procedure and allows fast recovery. It does not require general anesthesia so it can be done in your practice.

Minimal Downtime

Patients return home the same day and can resume normal activities within a few days.

The Sciton Edge

A Secure Investment

While laser-assisted lipolysis is rapidly evolving, ProLipo *PLUS* will not become obsolete. Sciton laser systems are enabled for tomorrow's technology and applications. The laser systems provide an easy upgrade path as procedures evolve.

Low Operating Costs

Sciton's non-proprietary fiber system protects you from unacceptably high per-case disposable costs. You can acquire fibers from a variety of sources – protecting you from vendors that charge high fees for their exclusive fibers.

Complete Aesthetic Solution

Other aesthetic treatment modules can be added to your laser, increasing your incremental revenue without the purchase of another full-priced system.

Practice Support

As part of our commitment to your success, we provide practice assistance through a ProLipo *PLUS* Practice Support Kit. The kit is a comprehensive set of materials that will enable you to easily incorporate ProLipo *PLUS* into your practice and to effectively market the ProLipo treatment to prospective patients.

Clinical Training

Your ProLipo *PLUS* purchase includes funds for training you on basic liposuction, as well as the proper use and operation of the laser and surgical technique for laser-assisted lipolysis.

Specifications

Wavelength	Nd:YAG 1064 nm and 1319 nm
Power (max)	40 W/1064 nm and 40 W/1319 nm
Repetition Rate (max)	50 Hz
Pulse Width	10 - 150 ms

Accessories

Fiber	600, 1000 micron (Open Source Fibers)
Cannula	14 GA: 200, 250, 300 mm 18 GA: 100, 150, 200 mm



visit us: www.ProLipoPLUS.com



925 Commercial Street, Palo Alto, CA 94303
 Phone: 888.646.6999 | www.sciton.com

©2009 Sciton, Inc. All rights reserved. ProLipo PLUS is a registered trademarks of Sciton Inc.

2600-029-00 Rev. A